

QuickSling - Arm sling

ART.NO 20118

QuickSling is a padded sling that is used for simple immobilisation of injured shoulders, elbows and wrists. The sling is easy to apply. Cut the desired length to suit each patient individually. The material is non-elastic and maintains its shape and width throughout the period of use. Simple to adjust using the Y-shaped Velcro fastening attached directly to the sling.

QuickSling can be applied in different ways to meet different needs. Available in black. QuickSling is supplied in a practical paper tray containing two rolls of 5.5 cm x 12 m sling strap and two bags of 14 Velcro straps. Sufficient for 14 double slings or 20 single slings.

Indications: For relieving shoulders/arms in the event of acute injury/trauma.

Contraindications: For patients with sinus carotid sensitivity avoid pressure from the orthosis in the carotid sinus area.

Material: 53% Polyurethane, 38% Polyester, 9% Nylon. Latex free.

Washing instructions: Machine wash, 40°. Close the Velcro fastening when washing. Ideally, use a washing net.



Applied as a double sling.



MD Class I CE

Size	Color	Art.no
Universal	Black	20118-00-50