

Sandra - Post surgical bra

ART.NO 31147

Firm and comfortable seamless post surgical bra. Good for use in connection with breast surgery as the knit, super-soft material protects delicate and sensitive skin. Also suitable as a soft prosthetic bra, with pockets for padding/breast form.

- Measured compression level equivalent to compression class 1.
- Opens at the back, with wide hook-and-eye, adjustable in three steps (4x3).
- Wide and adjustable shoulder straps at the back.
- Double material in the cup; pocket for padding/breast form.
- Soft, wide, seamless seam under the bust.

Indications: Post-surgical bra with a measured compression level equivalent to compression class 1. Intended for use after different types of breast surgery.

Contraindications: Infection. Arterial circulatory disorders. High risk of thrombosis.

Material: 79% Polyamide, 11% Elastane, 10% Polyester. Latex-free.

Washing instructions: Machine wash at 60°C. For optimal color durability, washing at 40°C is recommended. Can be tumble dried at low heat. Do not use bleach, chlorine bleach or fabric softener. Do not iron. Use of a laundry bag is recommended.

Size	Art.no Sandra white, small cup A-C	Art.no Sandra black, small cup A-C	Art.no Sandra white, big cup D-F	Art.no Sandra black, big cup D-F
S	31147-19-40	31147-10-40	31147-49-40	31147-40-40
M	31147-19-50	31147-10-50	31147-49-50	31147-40-50
L	31147-19-60	31147-10-60	31147-49-60	31147-40-60
XL	31147-19-70	31147-10-70	31147-49-70	31147-40-70
2XL	31147-19-80	31147-10-80	31147-49-80	31147-40-80

Size guide

cm	65	70	75	80	85	90	95	100	105	110	
UK	30	32	34	36	38	40	42	44	46	48	
Cup size	A										
	B	S (Cup A-C)		M (Cup A-C)		L (Cup A-C)		XL (Cup A-C)		2XL (Cup A-C)	
	C										
	D										
	E	S (Cup D-F)		M (Cup D-F)		L (Cup D-F)		XL (Cup D-F)		2XL (Cup D-F)	
	F										



Deep neckline, wide and adjustable closure, and adjustable shoulder straps in the rear.

