

# MAMMA MIA pelvic belt

ART.NO 50102

Mamma Mia is a pelvic belt that provides uniform, stable support across the pelvic joints. Mamma Mia also provides comfortable lift under the abdomen and is simple to adjust as the stomach expands. The increased width combined with the soft material makes it particularly pleasant to use. Colour: black.

- Combination of soft positioning belt and stable support strap
- Easy to adjust without affecting the placement
- Skin-friendly material: light and breathable

## Indications:

Model	Art.no	Height center front	Height center back
Mamma Mia	50102	16 cm	16 cm

Designed to counteract loosening of the joints/ symphysiolysis during or after pregnancy. Also suitable for non-pregnant people looking for stable sacroiliac support.

## Material:

Ventilated, breathable textile material that is gentle on the skin.  
Polyester 51%, Polyamide 17%, Rubber 12%, Polyurethane 8%, Cotton 7%, Elastane 5%.

## Washing instructions:

Machine wash, 40°. Close the Velcro fastenings. Use a washing net. Do not bleach, do not dry clean. Do not iron or tumble dry.

## Measurement:

Hip circumference in cm.



Size	Hip circumference in cm
Small	80-95
Medium	90-105
Large	100-115
X Large	110-125
XX Large	120-135

