

# Ballerina - Posture bandage

ART.NO 50103

Ballerina is a functional posture bandage that provides a comfortable reminder to maintain good posture in the upper back. Broad, soft straps across the shoulders cross over the back and fasten in front using Velcro. The broad, stable back section works effectively on the shoulder section without causing unpleasant pressure above the diaphragm. Ballerina also provides thoracic relief. The straps are simple to cut to the desired length. Flexible and comfortable brace that is ideal for wearing under everyday clothes. Color black.

**Indications:** Used for advanced thoracic kyphosis. Also suitable for collar bone fractures, following spinal stenosis operations, to achieve optimal breathing function, and in the case of compressed vertebrae.

**Material:** Ventilated, breathable textile material. Polyester 40%, Polyamide 37%, Polyethylene 17%, Elastane 5%, Elastomer 1%. Latex free.

**Washing instructions:** Machine wash, 40°. Close the Velcro. Use a washing net. Do not bleach, do not dry clean. Do not iron or tumble dry.

**Measurement:**

Thorax measurement below the bust in cm.

Size	Measurement in cm	Art.no
Small	65-80	50103-00-40
Medium	80-95	50103-00-50
Large	95-110	50103-00-60
X-Large	110-125	50103-00-70



MD Class I CE