

MammaFix - SI belt

ART.NO 50104

MammaFix is a semi-elastic SI belt (Sacroiliac Belt) that provides effective support across the pelvis. Suitable for use during pregnancy and in the case of pelvic joint problems. The belt is wrapped around the body and then fastened with Velcro in the centre. This makes it simple to adjust to provide the required stability. Sits securely in position as it is coated with silicone on the inside to prevent slipping. MammaFix is available in two lengths; both sizes are simple to shorten by cutting the belt at the marked seam along the edges of the strap. Made of sturdy elastic and rigid Velcro material. Colour: black.

Indications: Pregnancy-related pelvic girdle pain (PGP) or symphysis pubis dysfunction (SPD). Also suitable for non-pregnant people seeking stable support for the pelvic joints.

Contraindications: None known.

Material: 49% Polyester, 25% Natural Rubber Latex, 14% Silikon, 12% Polyamid.

Washing instructions: Machine wash, 40°C. Close the Velcro fastenings. Use a washing net. Do not bleach, do not dry clean. Do not iron or tumble dry.

Measurement: Hip circumference in cm.

Size	Hip circumference in cm	Art.no	Height
S/M	75-110	50104-00-40	7 cm
L/XL	110-145	50104-00-60	



LATEX Latex tread is covered by other textile fibers, will not get in contact with the skin.

