

AtleMage - Back orthosis

ART.NO 50105

AtleMage is a semi-elastic back brace that provides stable support across the lower back while simultaneously lifting/supporting the abdomen. Suitable for use by people with a large/loose abdominal section, as well as during and after pregnancy. The back section is made of AirX – an elastic material that provides good compression and is very comfortable to wear as it allows the skin to breathe. The front section is non-elastic and ensures stable lift/support for the abdomen. The 'shape' of the front section is adjustable on both sides, with Velcro fastenings allowing optimal lift/size. The brace is therefore also suitable for pregnant women and can be adjusted successively as their abdominal section expands. The elastic strap is secured in front and under the stomach with Velcro. It is adjustable to provide the required support across the lower back. Two replaceable back splints can be shaped individually to suit the contours of the back. Color beige.

Indications: Used for back pain, sciatica and other commonly occurring back problems. Ideal for providing extra support for a large/loose abdomen, during pregnancy and for people with minor hernias.

Material: Polyamide 46%, Polyester 44%, Elastane 5%, Cotton 4%, Natural Rubber Latex 2%.

Washing instructions: Machine wash, 40°. Close the Velcro fastenings. Use a washing net. Dry flat. Do not bleach, do not dry clean. Do not iron or tumble dry.

Measurement: Hip circumference in cm.

Size	Measurement Hip cm	Art.no	Height center front	Height center back
S	80-95	50105-08-40	25 cm	26,5 cm
M	95-105	50105-08-50		
L	105-115	50105-08-60		
XL	115-125	50105-08-70		
2XL	125-135	50105-08-80		



Latex tread is covered by other textile fibers, will not get in contact with the skin.

