

Tuff - Back orthosis

ART.NO 50200

Tuff back brace provides extremely stable support across the lower back/pelvis, and is ideal for people looking for a supple low back brace that provides strong support. The brace can be applied at waist height or over the pelvis, depending on where support is needed. Choose a smaller size if the brace is to be used at waist height. Tuff features broad Velcro straps on the front in order to ensure a snug fit over the hips or around the waist. The elastic straps are adjustable and provide stable support; they fasten front and centre with Velcro. Made of robust elastic and suits most figure shapes – for both women and men. Suitable for pain conditions in the lumbar spine, hip and pelvis. Also to ease residual pelvic pain after pregnancy. Tuff is available in black. Height center front 13-15 cm, height center back 20 cm.

Indications: Lumbar/lumbosacral back pain.

Contraindications: Unstable back fracture.
Sciatica pain of unclear genesis.

Material: Polyester 50%, Nylon 30%, Spandex 20%. Latex free.

Washing instructions: Machine wash at 40°C (104°F). Seal the Velcros. Use washing bag. Do not bleach. Do not dry clean. Do not iron. Do not tumble dry.

Measurement: Hip circumference in cm.

Size	Hip circumference in cm	Art.no
X-Small	65-80	50200-00-30
Small	80-95	50200-00-40
Medium	95-105	50200-00-50
Large	105-115	50200-00-60
X Large	115-125	50200-00-70
2X-Large	125-135	50200-00-80

