

Lyft - Binder

ART.NO 50206

Elastic abdominal binder which provides a comfortable and efficient support/compression of the abdomen. The front panel is shaped to "lift" the abdomen and give enough room to fit a pendulous figure. The elastic material is comfortable and feels soft on the skin. Color beige.

Level of compression: Medium, 2 out of 3.

Indications:

To be used post-operatively after abdominal surgery, to provide relief for surgical wounds and to prevent incisional hernia, also for pain relief and early mobilization of patients. To be used as an ostomy- or hernia support.

Counter indication: incarcerated hernia.

Material:

Polyester 77%, Natural Rubber Latex 17%, Polyamide 6%.

Washing Instructions:

Machine wash at 40°C (104°F). Seal the Velcros. Use washing bag. Do not bleach. Do not dryclean. Do not iron. Do not tumble dry.

Measurement:

To determine the correct size, measure hip circumference in centimeters.

Size	Measurement Hip cm	Art.no	Height center front	Height center back
S	80-95	50206-08-40	24 cm	23 cm
M	90-105	50206-08-50		
L	100-115	50206-08-60		
XL	110-125	50206-08-70		
2XL	120-135	50206-08-80		



Latex tread is covered by other textile fibers, will not get in contact with the skin.

