

StomiBas - Stoma binder

ART.NO 50210, 50211, 50212

Elastic stoma binder especially designed for stoma patients, to support the stoma and/or a parastomal hernia. The uniquely designed front panel is made of a non-elastic and non-fringing material and allows cutting a hole of any size and shape. The binder can also be used as traditional stoma belt without a hole. The product is available in three different heights, 16, 20 and 25 cm, in many sizes. The wide Velcro adjustable front closure makes the product easy to put on and adjust to requested compression. Can be closed right or left side.

Level of Compression: Strong, 3 out of 3.

Indications:

To be used as an ostomy support/hernia belt with or without hole for the stoma. The belt can also be used post-operatively after abdominal or thoracic surgery, to provide relief for surgical wounds and to prevent incisional hernia, also for pain relief and early mobilization of patients.

Material: Polyester 60%, Polyamide 27%, Rubber 13%.

Washing instructions: Machine wash at 40°C (104°F). Seal the Velcros. Use washing bag. Do not bleach. Do not dry clean. Do not iron. Do not tumble dry.

Measurement: To determine the correct size, measure circumference of the intended placement.

Size	Circumference cm	Art.no StomiBas 16cm	Art.no StomiBas 20cm	Art.no StomiBas 25cm
X Small	70-85	50211-08-30	50212-08-30	50210-08-30
Small	80-95	50211-08-40	50212-08-40	50210-08-40
Medium	90-105	50211-08-50	50212-08-50	50210-08-50
Large	100-115	50211-08-60	50212-08-60	50210-08-60
X Large	110-125	50211-08-70	50212-08-70	50210-08-70
2X Large	120-135	50211-08-80	50212-08-80	50210-08-80
3X Large*	130-145		50212-08-90	50210-08-90
4X Large*	140-155		50212-08-91	50210-08-91

*only available in 20 & 25 cm height.



50211: 16 cm



50212: 20 cm



50210: 25 cm

