

Bure - Back orthosis

ART.NO 50241, 50242

Stable back support in airy and thin material that provides extra good stability. Elastic straps and moldable stays for adjustable support. Suitable for acute backache, herniated disc or other pain conditions in the lower back. Suits most body shapes – men and women alike. Bure fastens in front with Velcro straps. The broad, elastic straps are adjustable, and they fasten front and centre with Velcro as well. The two replaceable back splints can be shaped to suit the contours of the back. Available in two heights.

Indications: Bure Low: Lumbar/lumbosacral back pain.

Bure High: Lumbal/lumbothoracal back pain.

Contraindications: Unstable back fracture.

Sciatica pain of unclear genesis.

Material: Polyester 55%, Polyamide 30%, Spandex 15%. Latex free.

Washing instructions: Machine wash at 40°C (104°F). Remove the splints. Seal the Velcros. Use washing bag. Do not bleach. Do not dry clean. Do not iron. Do not tumble dry.

Measurement: Hip circumference in cm.



50241: 15-17 cm front 26 cm back



50242: 19-21 cm front 32 cm back



MD Class I CE

Size	Measurement Hip cm	Bure Low, black			Bure High, black		
		Art.no	Height center front	Height center back	Art.no	Height center front	Height center back
Small	80-95	50241-00-40	15-17 cm	26 cm	50242-00-40	19-21 cm	32 cm
Medium	90-105	50241-00-50			50242-00-50		
Large	100-115	50241-00-60			50242-00-60		
X Large	115-125	50241-00-70			50242-00-70		
2X Large	125-135	50241-00-80			50242-00-80		