

Ronja - Back orthosis

ART.NO 50531, 50532, 50533, 50536

Stable fabric corset with elastic straps on the sides. Simple to put on and adjust for the required compression/stability – even for people with impaired hand strength. The fabric is slightly elastic, which ensures an excellent fit and makes the corset extremely comfortable to wear. The corset has two open elastic strips in the back and comes with two splints. Velcro fastening front and centre. Numerous thin plastic springs sewn in around the corset provide relief and stability. Color gray.

Ronja Semi-finished: Ronja Semi-finished – art.no 50536-26 – features two open seams in the back section and can be sewn in, let out and shortened as required. Edging, elastic strips and splints included.

Indications: Ronja low/Medium: Lumbar/lumbosacral back pain. Ronja High: Lumbal/lumbothoracal back pain.

Contraindications: Unstable back fracture. Sciatica pain of unclear genesis.

Material: 60% Cotton, 28% Polyamid, 7% Elastane, 5% Polyester. Latex free.

Washing instructions: Machine wash at 40°C (104°F). Remove the splints. Seal the Velcros. Use washing bag. Do not bleach. Do not dry clean. Do not iron. Do not tumble dry.

Measurement: Hip circumference in cm.



		Ronja Low - 50531			Ronja Mid - 50532			Ronja High - 50533			
	Size	Hip cm	Art.no	Height front	Height back	Art.no	Height front	Height back	Art.no	Height front	Height back
Child	XX Small*	60-75		18 cm	28 cm	50532-26-20	20 cm	33 cm		23 cm	38 cm
	X Small	70-85	50531-26-30			50532-26-30			50533-26-30		
	Small	80-95	50531-26-40			50532-26-40			50533-26-40		
	Medium	90-105	50531-26-50			50532-26-50			50533-26-50		
	Large	105-115	50531-26-60			50532-26-60			50533-26-60		
	X Large	115-125	50531-26-70			50532-26-70			50533-26-70		
	2X Large	125-135	50531-26-80			50532-26-80			50533-26-80		

*only model Ronja Mid - 50532

Ronja Semi-finished - 50536				
Size	Hip cm	Art.no	Height front	Height back
X Small-Small	70-95	50536-06-40	22 cm	37 cm
Medium-Large	95-115	50536-06-50		
X Large-2X Large	115-135	50536-06-60		

