

Frejka - Hipabduction orthosis

ART.NO 60232

The Hipabduction rail is used to position the hip joints of children having instability in the hip joint. Frejka consists of a flexible polyethylene rail with padded edges. Outside the rail pants of textile.

Indications: Congenital hip dislocation/subluxation. Hip dysplasia.

Material Composition: Polyethylene 49%, Cotton 29%, Polyamide 13%, Polyester 9%. Latex free.

Washing Instructions: Close all Velcro straps before washing. The Pant is machine washable at 40 degrees. Wipe the plastic part with a soft damp cloth, use mild detergent.



Size	Width, plastic part (cm)	Measurements between the knees (cm)	
X Small	14	15-18	
Small	16	17-20	
Medium	18	19-22	
Large	20	21-24	
X Large	22	23-26	
2X Large	24	25-28	



MD Class I CE