

# Knee-O-Trakker - Patella brace

ART.NO 70511

Knee-O-Trakker is a patella stabilizing orthosis with adjustable patella support. C-shaped neoprene pillow grips the skin and provides maximum patella control. The orthosis is put on with the pillow placed either medially or laterally. The pressure/support for patella is adjusted and controlled with two velcro straps.

**Indications:** Subluxation/luxation of patella and patellofemoral pain syndrome (PFPS).

**Contraindications:** None known.

**Size:** XS-4XL Universal right/left.

XS-S: children's sizes.


2XL-4XL: King-size.

**Material:** 50% Neoprene, 25% Nylon, 10% EVA Laminate, 5% Derlin, 5% Elastic, 5% Stretch UBL. Latex-free.

**Washing instruction:** Hand wash max 30°C. Dry flat. Secure all hook and loop closures before washing. Do not use bleach, chlorine bleach or fabric softener as this could impair the function of the material. Do not dry clean. Do not iron. Do not tumble dry.

**Measurement:** Circumference around knee, over patella.

	Size	Circumference cm	Art.no
Child	XS	30,5-33	70511-00-30
	S	33-35,5	70511-00-40
	M	35,5-38	70511-00-50
	L	38-40,5	70511-00-60
	XL	40,5-46	70511-00-70
King-size	2XL	46-51	70511-00-80
	3XL	51-56	70511-00-90
	4XL	56-61	70511-00-91



MD Class I CE



Fits both  
Adult & Child