

Dolly – Compression body

NEW



Adjustable shoulder straps

Adjustable hook and eye closure in 3 steps

Overlapping opening in the crotch

Shaped knitting over buttocks


NordiCare[®]

NORDICARE ORTOPEDI + REHAB AB • Solrosvägen 1 • SE-263 62 Viken, Sverige
Tel: 042-35 22 20 • Fax: 042-35 22 21 • info@nordicare.se • www.nordicare.se

Dolly – Compression body

NEW

Effective compression that provides fast and pleasant healing

- Recommended after surgery on buttocks, hips, abdomen and back.
- Shaped knitting over butt suitable after BBL and butt lift.
- Provides compression and support for pain relief and reduced swelling.
- Hook and eye closure in the middle front for easy application and adjustable fit.
- Knitted seamless material in 4-way stretch that shapes to the body.
- The hidden wrap opening in the branch provides compression in the area.



The adjustable shoulder straps can be attached in the middle front or sides.




Wide cut in the back and flexible fit.



Shaped knitting to follow the body contour.

Dolly, black

Size	Waist	Hip	Thigh	Art.no
XS	62-68 cm	86-91 cm	46-50 cm	51207-00-30
S	69-75 cm	92-97 cm	50-54 cm	51207-00-40
M	76-83 cm	98-103 cm	54-58 cm	51207-00-50
L	84-92 cm	104-111 cm	58-62 cm	51207-00-60
XL	93-101 cm	112-119 cm	62-66 cm	51207-00-70
2XL	102-110 cm	120-127 cm	66-70 cm	51207-00-80



Compression class 1



Adjustable with hook & eye



Seamless



4-way stretch

CE

Latex free

MD Class I

