

Next generation compression garments

A unique product assortment combining the most comfortable fabric together with high quality compression and functional design.



Seamless



Compression class 1



4-way stretch



40-60°C washing

Smart products, smart design

NordiCare's line of seamless knitted compression garments are designed to be worn directly after surgery and during the recovery period. To facilitate wear, we have incorporated a lot of smart design such as double closures, silicone waist line and open crotch. They are basically invisible under your clothing which helps you feel good while getting better. The line is constantly evolving and will in the near future expand with complementary models.

Seamless 4-way stretch



The knitted seamless 4-way stretch make the garments comfortable to wear and easy to fit as they adjust to body shape. The breathable and soft material transports moisture and is yet quick drying, durable and supportive. Small parts can be cut for drains or other facilities. The garments easily stay fresh and hygienic since they can be washed at home in high temperatures (40-60°C).

Targeted compression



Each garment is designed to provide targeted compression (class 1), meaning necessary compression and support exactly where needed. By controlling compression in the operated area, the healing process advances faster and patient comfort is heightened. For a perfect fit, the knitted zones and seams follow the contours of the body. Each product is specially designed based on indication and will ensure the right level of compression over time.



MD Class I




Jill

Compression bolero with full length sleeves. Recommended after arm and back surgery such as liposuction or arm lift. Hook and eye closure in the front for easy application and adjustment. Flexible bust area. Sizes XS-2XL. Art.no 51205.

Adjustable hook and eye closure in 3 steps

Fits different body shapes and bra sizes

Lighter compression over the abdomen



Alison



Compression body with covered back. Recommended after surgery or liposuction on the abdomen, sides and/or back. Supportive over the abdomen/pocket. Sizes XS-2XL. Art.no 51203.

Adjustable shoulder straps


Adjustable hook and eye closure in 3 steps both sides

Open crotch

CCL1
15-21 mmHg



Moveable shoulderstraps



High and wide back




Patsy

Full length compression pants. Recommended after surgery and/or liposuction on the leg or abdomen. Sizes XS-2XL. Art.no 51201.


Silicone waistband

Adjustable hook and eye closure in 3 steps both sides

Open crotch



CCL1
15-21 mmHg



Easy to apply




Annie

Below knee compression pant in pull-on model. Recommended after minor surgery or liposuction on the leg and/or abdomen. Sizes XS-2XL. Art.no 51202.

Silicone waistband

Can be used as step-2-product



CCL1
15-21 mmHg

Daisy

Fully openable compression body with double material over the abdomen/pocket. Recommended after surgery and/or liposuction on the abdomen. Sizes XS-2XL. Art.no 51204.

Silicone waistband



Adjustable hook and eye closure in 3 steps both sides

Crotch closure, hook and eye



Fully openable

CCL1
15-21 mmHg



Easy to apply



NordicCare®

NORDICARE ORTOPEDI + REHAB AB • Solrosvägen 1 • SE-263 62 Viken, Sverige
+46 (0)42-35 22 20 • info@nordicare.se • www.nordicare.se